decon
mobility

SMARTCHAIR FLOW

Dansk
English
Norsk
Svenska

Instructions for use

Indhold

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Contents

1. NOTICES	33
2. SAFETY INFORMATION	34
3. BATTERY SAFETY INFORMATION	35
4. SPECIFICATIONS	37
5. UNBOXING	39
6. UNFOLDING	40
7. FOLDING	41
8. SET UP	42
9. CONTROLLER	44
10. HOW-TO FEATURES	46
Backrest Height Adjustment	46
Armrest Height Adjustment	46
Joystick Forward and Backward Adjustment	46
Footrest Height Adjustment	46
Free Wheel	47
11. TRANSPORTING	48
12. BATTERY MAINTENANCE	50
13. CHARGING	52
14. MAINTENANCE, CARE, TROUBLESHOOTING	53
15. LIMITED WARRANTY	55
16. CONTACT & SUPPORT	56
17. APPENDIX 1 EMC DECLARATION	57
18. SYMBOLS	60

1. NOTICES

Thank you for choosing SmartChair Flow by Suzhou KD, a convenient and portable Powered Wheelchair.

It is important to read through this manual before operating your SmartChair Flow to become familiar with its features as well as its limitations and safety information.

You will find the operation and care of your SmartChair Flow very simple and easy, though should you have any questions, please contact your dealer or distributor if you need further assistance and guidance.

Indications for Use

It is a motor driven, indoor and outdoor transportation vehicle with the intended use to provide mobility to a disabled or elderly person limited to a seated position.

DISCLAIMER:

The notes, warnings, and cautions contained within this manual and marked by the triangular caution symbol should be given special attention. Users should also pay close attention to information marked in this manual with NOTICE.

Suzhou KD, its affiliates, and partners cannot be held responsible for personal injury or property damage resulting from unsafe or improper operation or maintenance of the SmartChair Flow power wheelchair. For questions about safe use, servicing, or maintenance, please contact your dealer or distributor.

As part of our ongoing product improvement initiative, Suzhou KD reserves the right to change specifications and design without notice. As a result, there may be minor differences between your SmartChair Flow and the photos, illustrations, or instructions in this manual.

2. SAFETY INFORMATION

- DO make certain that the seat is in a locked position when you are operating your SmartChair Flow.
- DO reduce your SmartChair Flow speed when you are driving around a corner.
- DO proceed with caution as you approach the downgrade of a ramp or other incline and when driving near a raised surface or unprotected ledges and drop offs (curbs, porches, stairs etc.).
- DO avoid exposure to rain snow, ice, salt, or standing water, whenever possible.
- DO maintain and store in a clean and dry condition.
- DO keep your hands away from the tires when driving. Be aware that loose fitting clothing can become caught in drive tires.
- DO keep tools and other metal objects away from battery terminals. Contact with tools can cause electrical shock.
- DO follow the recommended weight capacity (300lbs /136 kg).
- DO follow the recommended incline slope rating (10 degrees).
- DO check the battery charge to ensure there is enough charge to get to your destination before using.
- DO ensure your SmartChair Flow power is OFF when entering and exiting the power wheelchair.
- DO NOT drive your SmartChair Flow across the side of a hill or diagonally up or down a hill, and do not stop (if possible) while driving up or down an incline.
- DO NOT operate your new SmartChair Flow for the first time without completely reading and understanding the owner's manual.
- DO NOT carry additional passengers on your SmartChair Flow.
- DO NOT use your SmartChair Flow as a seat in a moving vehicle.
- DO NOT allow unsupervised children to play near the power wheelchair.
- DO NOT operate your SmartChair Flow while under the influence of alcohol, as this may impair your ability to drive safely.
- DO NOT sit in the SmartChair Flow while in a moving vehicle, as an accident or sudden stop may cause you to be thrown from the power wheelchair and compromise the seat lock latch.
- ALL servicing must be performed by qualified personnel only. Unauthorized modifications to power wheelchairs, its power system or its control system may compromise the performance and safety of the system. UNDER NO CIRCUMSTANCES should you attempt any repair or service that is not covered in this manual.
- ALWAYS ensure all safety warning labels are in place and legible.
- Persons with developmental impairments should not use the device unsupervised, and users should not use the device while impaired in any way.
- The device is for passenger use only. Please do not allow children to use it.
- Please place the device in a well-ventilated area that is easy to access while charging.
- If the wheelchair loses control, please immediately disconnect the battery or turn off the controller's power to prevent any potential harm.
- The wheelchair cables are not toys. Please do not wrap them around your neck to avoid the risk of suffocation.

⚠WARNING: The operator should read all the warning signs and instructions here and on the power wheelchair before use.

⚠WARNING: You should not operate your power wheelchair on public streets and roadways. Be aware that it may be difficult for traffic to see you when you are seated on your power wheelchair. Obey all local pedestrian traffic rules. Wait until your path is clear of traffic and then proceed SmartChair Flow with extreme caution.

3. BATTERY SAFETY INFORMATION

WHEN USING THE BATTERY

1. Lithium-Ion battery packs may get hot, explode or ignite and cause serious injury if exposed to abuse conditions. Be sure to follow the safety warnings listed below:
 - DO NOT place the battery near fire or heat the battery.
 - DO NOT install the battery backwards so the polarity is reversed.
 - DO NOT connect the positive and negative terminals of the battery to each other with any metal object(such as a wire).
 - DO NOT carry or store the battery together with necklaces, hairpins, or other metal objects.
 - DO NOT pierce the battery with nails, strike the battery with a hammer, step on the battery or otherwise subject it to strong impacts or shocks.
 - DO NOT solder directly on to the battery.
 - DO NOT expose battery to water or salt water, or allow the battery to get wet.
2. DO NOT disassemble or modify the battery. The battery contains safety and protection devices, which if damaged, may cause the battery to generate heat, explode or ignite
3. DO NOT place the battery in or near fire, on stoves or other high temperature locations. Do not place the battery in direct sunlight, or use or store batteries inside cars in hot weather. Doing so may cause the battery to generate heat, explode or ignite. Using the battery in this manner may also result in a loss of performance and a shortened life expectancy.
4. The device is for adult use only.
5. When the battery is worn out, insulate the terminals with adhesive tape or similar materials before disposal.
6. Immediately discontinue use of the battery if, while using, charging or storing the battery, the battery emits an unusual smell, feels hot, changes color, or shape, or appears abnormal in any way. Contact your dealer or distributor if any of these problems are observed.
7. DO NOT place battery in microwave oven, high pressure containers or on induction cookware.



CHARGING PRECAUTIONS

1. BE SURE to follow the warnings listed below while charging the battery. Failure to do so may cause the battery to become hot, explode or ignite and cause serious injury.
 - When charging the battery, ONLY USE THE SPECIFIED CHARGER PROVIDED WITH THE PRODUCT.
 - DO NOT attach the battery to a power supply plug or directly to a car's cigarette lighter.
 - DO NOT place the battery in or near fire, or in direct sunlight. Heating the battery can damage the safety circuitry, which can cause additional heating, rupture or ignition of the battery.
2. DO NOT continue charging the battery if it does not recharge in the timelines specified in the operating manual (6-7 hours or until the charging brick turns green).
3. THE TEMPERATURE range over which a battery can be charged is 0°C to 45°C (32°F to 113°F). Charging the battery at temperatures outside of that range may cause severe damage to the battery or reduce battery life expectancy.

WHEN DISCHARGING THE BATTERY

1. DO NOT discharge the battery using any other device except for the specified device. When the battery is used in devices other than the specified device, it may damage the battery or reduce its life expectancy. The device may cause abnormal current flow, which can cause the battery to become hot, explode or ignite and cause serious injury.
2. THE TEMPERATURE range over which the battery can be discharged is -5°C to 45°C. Use of the battery outside this temperature range may damage the performance of the battery or may reduce its life expectancy.

4. SPECIFICATIONS


DIMENSIONAL SPECIFICATIONS

UNFOLDED DIMENSIONS:	980mm*600mm*908mm
FOLDED DIMENSIONS:	750mm*600mm*350mm
CHAIR GROSS WEIGHT:	25KG
CHAIR NET WEIGHT:	23KG (NO BATTERY)
SEAT CUSHION SIZE:	460mm(W)*450mm(D)
SEATING WIDTH:	460mm
SEATING DEPTH:	470mm
FRONT WHEEL DIAMETER:	8"
REAR WHEEL DIAMETER:	10"
SEAT SURFACE HEIGHT AT FRONT EDGE:	480mm
BACK SUPPORT HEIGHT:	500mm/550mm/600mm
SEAT PLANE ANGEL	6°
BACKREST ANGEL	11°
FOOTREST TO SEAT DISTANCE	470mm
LEG TO SEAT SURFACE ANGEL	95°
ARMREST TO SEAT DISTANCE	200mm-300mm
FRONT LOCATION OF ARMREST STRUCTURE	340mm

PERFORMANCE SPECIFICATIONS

DRIVING RANGE (energy consumption):	25 km
MAXIMUM SPEED:	5.868 km/h (1.63 m/s) Q3
DRIVE SYSTEM:	Direct drive on the rear wheels
MOTOR TYPE:	Brushless DC24V250W, 2pcs
INCLINE CAP ABILITY (Dynamic stability uphill):	10°
TURNING RADIUS:	800mm
BRAKING DISTANCE:	≤1 m
MAX OBSTACLE CLIMBING:	50mm
BRAKING SYSTEM:	ELECTROMAGNETIC
MAX WEIGHT CAPACITY:	136KG
STATIC STABILITY DOWNHILL	Min 25.7° max 27°
STATIC STABILITY UPHILL	Min 20.4° max 22.2°
STATIC STABILITY SIDEWAYS	Min 19.3° max 19.5°
BATTERY TYPE	Lithium-Ion (24VX12AH)
BATTERIES PROVIDED:	ONE (1)
BATTERY WATT HOUR:	288WH
HORN:	BUILT-IN
TIRES:	PU
CHARGING TIME:	6 TO 7 HOURS
CHARGER INPUT:	100-240V,50/60HZ(INTERNATIONAL)
CHARGER OUTPUT:	24V 2.0A
LAP BELT	Provided

5. UNBOXING


- To remove the powerchair, orientate the box, so it is laid down on the ground.
- Grab the top and begin to pull your SmartChair Flow out horizontally parallel with the ground.

- Locate the components provided with your power wheelchair

Seat cushion

Lithium battery


Joystick

Battery charger


Tool kit

Tool

6. UNFOLDING

1. Lift the footrest and pull the safety lock upward to unlock the chair.

2. Unfold the wheelchair outward until you hear a clicking sound indicating it is locked in place.

3. Lower the anti-tip wheel bracket, complete the unlocking of the wheelchair.

⚠ NOTICE: Failure to fully extend the power chair into its locked position may result in injury while driving.


7. FOLDING

1. Lift both the armrests and the anti-tip wheel brackets upward.


2. Pull the unlocking handle and push the wheelchair forward.


3. Until you hear a clicking sound indicating it is locked in place.

8. SET UP



- Slide the joystick controller into the armrest tube until a comfortable arm length is achieved for the user (image 1).
- Align the grey points of the joystick connection cables and press them together to plug the joystick controller in (image 2).
- Rotate the coupling over the connection point to secure the connection point.
- Make sure the quick-release pin locks into the tube to complete the fixing action.

- Before sitting in the electric wheelchair, be sure the power button is OFF and the rear brake levers are in the locked position.
- SmartChair Flow offers the user 2 optional entry alternatives, direct front entry or side entry
- For a front entry, we recommend flipping up the foot plate to create an easier sitting experience.
- DO NOT step on the footplate when entering or exiting the SmartChair Flow as it may cause the power chair to move in an unsafe manner.

- To enter on the side, lift the arm rest vertical to provide the required clearance.
- To enter on the side, lift the arm rest vertical to provide the required clearance.

- Joysticks can be positioned on left or right arm rest.
- Rotate the joystick connector to loosen it (Image 1).
- Pull up the quick-release pin to remove the joystick (Image 2).
- Insert the joystick fixing plate into the armrest tube on the other side. Adjust to the desired position and lock the quick-release pin (Image 3).
- Insert the joystick connector (Image 4) and rotate the locking ring to secure it.
- Ensure the quick-release pin is fully locked into the groove to complete the installation.
- Move the cable to the new side and secure it with ties. Make sure it doesn't block folding or movement, then reconnect the joystick.

9. CONTROLLER

- 1. Power On/Off Button** – Press this button briefly to power on the powered wheelchair and hold it down to power off.
- 2. Speed Display** – Displays the current driving speed of the wheelchair in real time. As you accelerate or decelerate, the numbers on the screen will change accordingly, helping you easily monitor your current speed.
- 3. Speed Indicator**—One bar represents the slowest speed, and five bars represent the fastest speed.
- 4. Charging Indicator**— When the device is connected to the charger, the light will blink during charging; The indicator light will remain steadily on after charging is complete.
- 5. Battery Display** – Displays the real-time battery level of the powered wheelchair.
- 6. Mode Display**— There are a total of three modes: SOFT, NORMAL, SPORT.
- 7. Decelerate Button** – Press this button to decrease speed. Hold the brake button down, and the driving mode of the wheelchair will cycle from SPORT to NORMAL, and finally to SOFT. You can adjust based on your preference
- 8. Accelerate Button** – Press this button to increase speed, with a maximum speed of 6 kilometers per hour. Hold the accelerate button down, and the driving mode will change from SOFT to NORMAL, and finally to SPORT.
- 9. Horn Button** – Press this button to activate a beeping sound from the powered wheelchair to alert the surrounding environment.
- 10. Joystick** – Controls the direction and speed of the powered wheelchair; push the joystick in the direction you wish to go.

11. **Type-C Power Supply Port** – Use a power cable with a 5V/2A output specification to connect mobile devices for charging.
12. **Wheelchair Charging Port**– You can charge the powered wheelchair here without removing the battery.

10. HOW-TO FEATURES


1.

2.

Backrest Height Adjustment


By loosening the four set screws on the top of the backrest (Image 1), the height of the backrest can be freely adjusted. The adjustable range is 100mm (Image 2).


3.

Armrest Height Adjustment

By loosening the quick-release handle on one side of the armrest, you can freely adjust the height of the armrest easily, the adjustable range is 100mm. based on your preference (Image 3).

4.

5.

Joystick Forward and Backward Adjustment


Using the quick-release button located at the base of the armrest (Image 4), you can easily pull it to quickly adjust the armrest to three different positions (Image 5). The total adjustable range is 50mm.

6.




Footrest Height Adjustment

By loosening the two screws located at the footrest (Image 6), you can adjust the height of the footrest to three different levels. The adjustable range is 50mm.

Push mode

- In order to have a care giver push the SmartChair Flow, it must be put in PUSH mode
- This mode disengages the drive mechanism to the wheels, making them free to push manually.

Lower anti-tip wheels

- Please ensure the anti-tip wheel brackets are in the lower position during use to help reduce the risk of tipping and improve overall safety.
- Press the anti-tip wheel bracket downward; the spring pin will automatically slide into the built-in round recessed hole on the motor (Image 1).
- Securely locking the anti-tip wheel into position (Image 2).

Fold anti-tip wheels


- Pull the anti-tip wheel unlock button (Image 3). The anti-tip wheel will automatically fold upward via the internal tension spring.
- It will then automatically lock into place through the spring mechanism within the unlock button (Image 4&5)

11. TRANSPORTING

Towing Function Guide


After folding the wheelchair, you can easily tow your powered wheelchair by following these steps:

- Flip up the footrest (Image 1).
- Grab the hollow part of the footrest (Image 2), and using the rear wheels, you can easily move your powered wheelchair.


Loading into the Trunk Function Guide

After folding the wheelchair, you can easily load it into the trunk of your car by following these steps:



1. Move the powered wheelchair to the rear of your car's trunk (Image 1), ensuring that the Smartrolls are positioned close to the trunk.
2. For user convenience, the Smartrolls are designed with two adjustable positions to accommodate both sedans and full-size SUVs. You can select the appropriate angle based on the height of your vehicle's trunk.

3. Pull the unlock pin to open the Smartrolls (Image 2, Image 3). The unlock pin has two different angles to accommodate different trunk heights. After adjusting to the first angle, you can continue adjusting to a larger angle.

4. Grasp both sides of the bottom of the backrest tube (Figure 8.2.4) so that the Smartrolls contact the flat surface of your car. Then, hold both sides of the bottom of the backrest frame (Figure 8.2.5) and lift your powered wheelchair upwards. Push it towards the interior of your car's trunk to successfully load it into the trunk.

12. BATTERY MAINTENANCE

Lithium-Ion rechargeable batteries require routine maintenance and care in their use and handling. Read and follow the guidelines in this document to safely use Lithium-Ion batteries and achieve the maximum battery life span.

OVERVIEW

Do not leave batteries unused for extended periods of time, either in the product or in storage. When a battery has been unused for 6 months, check the charge status and charge or dispose of the battery as appropriate.

Rechargeable Lithium-Ion batteries have a limited life and will gradually lose their capacity to hold a charge. This loss of capacity (aging) is irreversible. As the battery loses capacity, the length of time it will power the product (run time) decreases.

Lithium-Ion batteries continue to slowly discharge (self-discharge) when not in use or while in storage. Routinely check the battery's charge status. Use only the supplied battery in your product.

Observe and note the time that a new fully charged battery provides for powering your product. Use this new battery run time as a basis to compare run times for older batteries. The running time of your battery will vary depending on the product's configuration and the applications that you run.

Routinely check the battery's charge status. Carefully monitor batteries that are approaching the end of their estimated life.

Consider replacing the battery with a new one if you note either of the following conditions: (The battery can be replaced only by certified repair shops with trained professionals. Contact the authorized local Dealer or Distributor for assistance, alternatively direct contact with the manufacture.)

- The battery run time drops below about 80% of the original run time.
- The battery charge time increases significantly.

If a battery is stored or otherwise unused for an extended period, be sure to follow the storage instructions in this document. Not follow the instructions, and the battery has no charge remaining when you check it, consider it to be damaged. Do not attempt to recharge it or to use it. Replace it with a new battery.

CHARGING

Always follow the charging instructions provided with your product. Refer to your product's user manual and/or online help for detailed information about charging its battery.

Use only the supplied charger to charge your battery.

STORAGE

- Charge or discharge the battery to approximately 50% of capacity before storage.
- Charge the battery to approximately 50% capacity at least once every six months.
- Remove the battery and store it separately from the product.
- Store the battery at temperatures between 5 °C and 20 °C (41 °F and 68 °F).

HANDLING PRECAUTIONS

- Do not disassemble, crush, or puncture a battery.
- Do not short the external contacts on a battery.
- Do not expose a battery to fire or water.
- Do not expose a battery to temperatures above 60°C (140 °F).
- Keep the battery away from children.
- Avoid exposing the battery to excessive shock or vibration.
- Do not use a damaged battery.
- If a battery pack has leaking fluids, do not touch any fluids.
- Dispose of a leaking battery pack (see Disposal and Recycling in this document)
- In case of eye contact with fluid, do not rub eyes. Immediately flush eyes thoroughly with water for at least 15 minutes, lifting upper and lower lids, until no evidence of the fluid remains. Seek medical attention.
- Immediately discontinue use of the battery if, while using, charging or storing the battery, the battery emits an unusual smell, feels hot, changes color, or shape, or appears abnormal in any way. Contact your dealer or distributor if any of these problems are observed.

TRANSPORTATION

- Always check all applicable local, national and international regulations before transporting a Lithium-Ion battery.
- Transporting an end-of-life, damaged, or recalled battery may, in certain cases, be specifically limited or prohibited.

DISPOSAL AND RECYCLING

Lithium-Ion batteries are subject to disposal and recycling regulations that vary by country and region. Always check and follow your applicable regulations before disposing of any battery. Lithium-ion batteries are subject to disposal and recycling regulations that vary by country and region. In the EU, batteries must be disposed of through approved collection systems in accordance with the WEEE directive. Please contact your local battery recycling authority for guidance.

Many countries prohibit the disposal of waste electronic equipment in standard waste receptacles.

Place only discharged batteries in a battery collection container. Use electrical tape or other approved covering over the battery connection points to prevent short circuits.

13. CHARGING

Battery Charging Methods

SmartChair Flow is equipped with dual charging ports, located beneath the controller (Image 1) and at the bottom of the battery (Image 2). You can choose either charging method according to your needs.

1.

2.

- Fully Charged
- Charging

- When charging, the indicator light on the power charger is red.
- After 6 to 7 hours of charging (depending on the remaining battery level), the indicator light will change from red to green, indicating that the battery is fully charged.
- To check your battery level, turn ON the joystick controller power.
- There are a total of five WHITE battery indicator lights.
- From top to bottom, the lights will gradually turn off as the battery drains, indicating the remaining charge.
- Recharge the battery before the last WHIT light turns off, as this indicates the battery is fully depleted.

BATTERY CHARGE INDICATOR ON JOYSTICK CONTROLLER

14. MAINTENANCE, CARE, TROUBLESHOOTING

Joystick Malfunction Code

Code	Information
E1	Upper and lower control communication failure
E2	Push handle to power on
E3	Overvoltage
E4	Upper control joystick abnormality
L1	Left motor electromagnet abnormality
L2	Left motor phase line abnormality
L3	Left motor Hall sensor abnormality
L4	Left motor phase line short circuit
J1	Right motor electromagnet abnormality
J2	Right motor phase line abnormality
J3	Right motor Hall sensor abnormality
J4	Right motor phase line short circuit

It is important to perform regular maintenance of your Smartchair Flow for longevity and warranty eligibility. Even a little cleaning now and again will go a long way to keeping your device looking and functioning like new.

Certified repair shops with trained professionals who specialize in electric wheelchairs will be qualified to perform the service and maintenance on your Smartchair Flow electric wheelchair should it be required. Please be sure to keep a log of the service and maintenance performed. Your shop or technician may require parts, which your dealer or distributor has available.

CLEAN WHEELCHAIR

Please clean the wheelchair regularly. Clean the parts which are frequently touched with user's body (such as seat cushion, armrest, controller) by a clean and slight wet cloth. Do not use organic solvents to clean.

We recommend you inspect your SmartChair Flow two to three times per year for in the following areas:

Frame: Tighten loose screws as necessary using the tool kit included with your Smartchair Flow

Wheels: Make sure they rotate easily

Significant resistance may require new bearings

Armrests: Look for excessive wear or if they are loose

Replacement parts and more major maintenance

issues should be done by a professional. We do

not recommend that you do intense types of

maintenance work for the Smartchair Flow yourself.

LUBRICATION

- According to use condition, please clean or lubricate the front axle, front fork and rotating parts after 180 days
- Molybdenum disulfide grease is recommended
- Your Smartchair Flow Powered Wheelchair has electric hub motors that do not require lubrication, please do not lubricate the motors

Question	Possible Causes	Answer
Why is there NO POWER to the SmartChair Flow?	<ul style="list-style-type: none"> • Controller system power is not connected. • Controller system CPU and joystick are not connected. • Battery power is too low. 	<ul style="list-style-type: none"> • Connect the battery. • Ensure all connections between joystick, CPU and battery are securely tightened. • Charge the battery.
Why does my power chair beep and not move?	<ul style="list-style-type: none"> • Freewheel mode engaged. • Capacity is reached. • Other error code. 	<ul style="list-style-type: none"> • Ensure the rear motor switches are in POWER mode, rear wheels should be locked. • Contact your dealer or distributor if other error code present.

Working and storage temperature, humidity, and atmospheric pressure.

Working Environment Humidity: 20%-80%

Working Temperature: -5C°-45C°

Storage Environment Humidity: 20%-80%

Storage Temperature: - 10C°~ 40C°

Atmospheric Pressure: 860hpa~1060hpa

It takes 6 hours for the wheelchair to return to a usable condition after being exposed to the minimum or maximum temperature. Please allow it to reach room temperature before use.

Please do not place the device in environments with cotton, dust, direct light, or other factors that may affect the device, to avoid causing damage.

The expected service life of the device is 5 years.

15. LIMITED WARRANTY

The purchaser of this product is entitled to a limited warranty as offered by Suzhou KD and its affiliates on the following components and timeframes from the date of purchase:

5 YEAR WARRANTY - Main Frame

1 YEAR WARRANTY - Motor

1 YEAR WARRANTY -Control System

1 YEAR WARRANTY -Battery

90 DAY LIMITED WARRANTY

On wear parts, such as the tires, armrest padding, joystick button covers, joystick knobs, seat cushion, and backrest padding.

This warranty is valid for the replacement of dysfunctional parts only. Any parts under warranty will be replaced and shipped to your door. Any service and labor fees, if applicable, to replace parts under warranty must be paid by the user.

Due to its straightforward design, most parts can be easily exchanged by the end user without a professional service tech required. However, it is always recommended you seek professional help with maintenance and service, to make sure the work is done properly.

The warranty only covers NON-DELIBERATE Damage.

The warranty does not cover:

- Products damaged by user negligence.
- Products were damaged accidentally.
- Products damaged intentionally.
- Products that have been subjected to negligence.
- Products that have been subjected to abuse.
- Products that have been improperly stored.
- Products that have been improperly handled.
- Products that have been improperly operated.
- Products that have experienced general misuse.
- Products that have been modified in an unauthorized, unapproved way.

The warranty is non-transferable and only valid for the original wheelchair purchaser.

Suzhou KD will, at its sole discretion, replace any warrantable item that fails within the warranty period.

Contact your dealer and distributor if any warrantable item that fails within the warranty period

16. CONTACT & SUPPORT

You should report any suspected serious incident that has occurred in relation to the SmartChair Flow to the manufacturer, its authorised representative, the importer and the competent authority of the Member State in which you are established.

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17. APPENDIX 1 EMC DECLARATION

Warning: Don't use the device near active high frequency (HF) surgical equipment and the radiofrequency (RF) shielded room of an medical equipment (ME) system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Warning: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation."


Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the equipment, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Guidance and manufacturer's declaration - electromagnetic emissions.		
The device is intended for use in the electromagnetic environment specified below, The customer or the user of the device should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	Note: The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	CLASS B	Note: The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class B). Pst and Plt requirements shall not be applied to voltage changes caused by manual switching. The limits shall not be applied to voltage changes associated with emergency switching or emergency interruptions. Please also refer to Annex A (Application of limits and type test conditions) for details in IEC/EN61000-3-3.
Harmonic emissions IEC 61000-3-2	CLASS A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Clause 5	

Guidance and manufacturer's declaration - electromagnetic Immunity		
Immunity Test	IEC 60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±2 kV, ±4 kV, ±6 kV, ±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air ±2 kV, ±4 kV, ±6 kV, ±8 kV VCP, HCP	±2 kV, ±4 kV, ±6 kV, ±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air ±2 kV, ±4 kV, ±6 kV, ±8 kV VCP, HCP
Electrical fast transient/burst IEC 61000-4-4	±1 kV, ±2 kV for power supply lines ±1 kV signal input/output 100 kHz repetition frequency	±2 kV for power supply lines Not Applicable 100 kHz repetition frequency
Surge IEC 61000-4-5	±0.5 kV, ±1 kV live to neutral	±0.5 kV, ±1 kV live to neutral
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT for 0.5T At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT for 1T At 0° 0 % UT for 300T At 0°, 180°. (100V/60Hz) 0 % UT for 250T At 0°, 180°. (240V/50Hz) 70 % UT for 30T At 0° (100V/60Hz) 70 % UT for 25T At 0° (240V/50Hz)	0 % UT for 0.5T At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT for 1T At 0° 0 % UT for 300T At 0°, 180°. (100V/60Hz) 0 % UT for 250T At 0°, 180°. (240V/50Hz) 70 % UT for 30T At 0° (100V/60Hz) 70 % UT for 25T At 0° (240V/50Hz)
Proximity Magnetic Fields IEC 61000-4-39	8 A/m (for 30kHz) 65 A/m (for 134.2kHz) 7.5 A/m (for 13.56 kHz)	8 A/m (for 30kHz) 65 A/m (for 134.2kHz) 7.5 A/m (for 13.56 kHz)
Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC61000-4-6	3 V 0,15 MHz – 80 MHz	3 V 0,15 MHz – 80 MHz

Radiated RF IEC61000-4-3	6 V in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	6 V in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz
	20 V/m 26 MHz – 2.5 GHz (for Forward mode, Backward mode, Brake mode & Standby mode); 10 V/m 80 MHz – 2.7 GHz (for Wheelchair charging mode-main body charging port, Battery-pack charging mode, Forward mode, Backward mode, Brake mode & Standby mode) 80 % AM at 1 kHz	20 V/m 26 MHz – 2.5 GHz (for Forward mode, Backward mode, Brake mode & Standby mode); 10 V/m 80 MHz – 2.7 GHz (for Wheelchair charging mode-main body charging port, Battery-pack charging mode, Forward mode, Backward mode, Brake mode & Standby mode) 80 % AM at 1 kHz
NOTE UT is the a.c. mains voltage prior to application of the test level.		

18. SYMBOLS

<p>IP24</p>	<p>waterproof and dustproof resist</p>
	<p>Body Floating Equipment</p>
	<p>Limit of Humidity</p>
	<p>Waste Electrical and Electronic Equipment</p>
	<p>Medical Device</p>
	<p>Consult Instructions for Use</p>
	<p>Do not use the device after it is damaged</p>
	<p>General Warning Sign</p>
	<p>Manufacturer</p>
	<p>Unique Device Information</p>

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Decon Wheel AB