### INTRODUCING EPPUR DREEFT DISTRIBUTED BY



#### Sales contacts

Germany info@deconmobility.de +49 (0) 3222 18 52 896 Netherlands info@mobilityproducts.nl +31 (0)85 - 87 86 144 Scandinavia sales@decon.se +46 (0) 34540880



## \_DREEFT

DO YOU WANT TO BRAKE FREE?



Ξ P P U R

The first braking device for manual wheelwhairs


A SET OF WHEELS TO SLOW DOWN WITH NO PAIN

b co

Today, braking or turning a wheelchair implies using your hands as brake pads, causing shoulder pain and hand burns. Thanks to Dreeft wheels, brake with your fingertips with no effort or friction!

Just lightly pull the handrim backwards to activate the coaster brake system.





ECHNICAL DATASH	IEET	Spokes:	24 spokes 14G black steel
Vheel dimensions :	24" (ETRTO 540mm)	Handrim fixation part :	Black powder coated aluminium
yres:	Schwalbe Marathon+ 25-540	Quick release axis diameter :	Compatible with 12mm & 12,7mm
landrims :	Black anodized aluminium Ø19mm	Handrim-to-rim spacing :	12mm or 19mm
tims :	Double wall, aluminium	Max user weight:	120 kg
landrim grip :	PVC rubber hand grip Ø30mm	Fully equipped wheel weight:	: 3 750 g



#### Zero friction

Say goodbye to gloves and hand burns



#### **Motorizable**

Helps you steer with rear power assist. No pain, no danger!



#### 80% Less effort

Compared to standard wheels (Your shoulders will say thank you)



# Appropriate for para & quadri

Usable without hand gripping, even under the rain!