







**Meyra - Avanti Pro 1.735, Avanti 1.736, Basic 1.751
MED3052**







1	MED2000	1	
2	100 400 017	1	
3	100 400 007	1	
4	HW069	1	
5	XSM16M	1	
6	XSM8ML	1	
7	XSM8X45	1	
8	100 400 005	2	
9	100 400 012	2	
10	405242	1	




Code 81




MED1428L/R



1	MED1428R	1	
---	----------	---	---

 = MED1428L


Code 21, 84 & 748
MED4054R



1	XSM6BR	2	
2	XSM6ML	2	
3	XSM6X14	2	
4	MED1500	1	 
5	XOC-48160-B0	1	

 = MED4054L


Code 21, 84 & 748
MED4055R



1	XSM6BR	2	
2	XSM6ML	2	
3	XSM6X14	2	
4	X0F-48160-D00	1	
5	MED1500	1	

 = MED4055L


Code 107
MED4053R



1	MED1590R	1	 
2	XSM6X16	2	
3	XOC-48160-B0	1	
4	XSM6ML	2	

 = MED4053L

MED1218




1	XSM6BRTG	1	
2	XOD-2612K-00	1	
3	XA2-2612E-00	1	
4	XSM6BR	1	
5	MC6SM6X50	1	



MED1218 + MED1643



M3	1,2 Nm
M4	3 Nm
M5	6 Nm
M6	10 Nm
M8	25 Nm
M10	50 Nm
M12	85 Nm



M16x1	30 Nm
M18x1	40 Nm