

















**Ottobock - Start M2S Flex, Start M2S Flex Prop, Start M2S Flex Strong  
MED3330**












1	HW084	1	
2	MED2000	1	
3	405241	1	
4	XSM8X40	1	
5	100 400 014	1	
6	XSM18M	1	
7	100 400 008	2	
8	100 400 013	1	
9	100 400 015	1	
10	XSM8ML	1	

















1	MED4046	1	
2	XSM6ML	2	
3	XSM6BR	2	
4	XSM6X14	2	
5	XOC-48160-B0	1	
6	XSM6X35	2	

## MED1670




1	MED4046	1	
2	XSM6ML	2	
3	XSM6BR	2	
4	XSM6X14	2	
5	X0F-48160-D00	1	
6	XSM6X35	2	





1	MED1671	1	
2	MED1672	1	
3	XSM6BRTG	1	
4	XSM6BR	1	
5	MC6SM6X50	1	




=MED1673+MED1643



1	100 400 039	1	
2	XSM6X25	1	



M3	1,2 Nm
M4	3 Nm
M5	6 Nm
M6	10 Nm
M8	25 Nm
M10	50 Nm
M12	85 Nm



M16x1	30 Nm
M18x1	40 Nm