







**Invacare Clematis Pro
MED3359**







1	MED2000	1	
2	100 400 004	1	
3	100 400 005	1	
4	100 400 017	1	
5	100 400 027	1	
6	100 400 020	1	
7	XSM6X35	1	
8	405239	1	
9	HW069	1	
10	XSM16M	1	
11	XSM6ML	1	






MED1414L or MED1414R




1	MED1661 / MED1662	1	
2	XOC-48160-B0	1	
3	XSM6X16	2	
4	XSM6ML	2	
5	XSM6BR	2	
6	XSM4X20	2	
7	XSM4BR	2	



1	XSM6BRTG	1	
2	Xklamma2	1	
3	XSM6BR	1	
4	XSM6X20	1	




MED1252 + MED1643



1	XA2-2612K-00	1	
2	XA2-2612E-00	1	
3	XSM6BRTG	1	
4	MED1317	1	
5	XSM6BR	3	
6	MC6SM6X35	1	
7	XSM6X16	2	
8	MED1316	1	
9	XSM6ML	2	



M3	1,2 Nm
M4	3 Nm
M5	6 Nm
M6	10 Nm
M8	25 Nm
M10	50 Nm
M12	85 Nm



M16x1	30 Nm
M18x1	40 Nm