


Assembly instruction for Easywheel, Etac Cross 6 VH3010-75

1. Remove the original wheels and hub sleeves on both sides.
2. Install the new hub sleeves fitting Easywheel into existing parts.



1. Details for mounting on one side.

<i>Item Number</i>	<i>Part number</i>	<i>Qty</i>	<i>Description</i>
1	VH51001-04	1	Hub sleeve
2	XSM16M	1	Nut M16x1
3	405242	1	Lock-washer M16

2. All positions in the illustration below are available.

