





| Item No. | Part Number | Qty | Description         |
|----------|-------------|-----|---------------------|
| 1        | VH51002-03  | 1   | Navhylsa M18        |
| 2        | 405240      | 1   | Nordlock Bricka M18 |
| 3        | VH53001-02  | 2   | Distans             |
| 4        | VH52001-01  | 1   | Lock Platta         |
| 5        | XSM6X30     | 1   | SKruv M6 x 30       |
| 6        | XSM6ML      | 1   | Låsmutter M6        |
| 7        | XSM18FJ     | 1   | Fjäder Bricka M18   |
| 8        | XSM18M      | 1   | Mutter M18 x1       |

|                                                                                                                                |                                                                                                                                                                         |                                                    |                                                                                     |                      |                                                       |                                |
|--------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------|-------------------------------------------------------------------------------------|----------------------|-------------------------------------------------------|--------------------------------|
| <div>A4</div> <div></div> <div>ISO E</div> | <div>IronCAD</div> <div>2D-ritning</div> <div>3D-fil tillgänglig</div> <div></div> | <div>Generella</div> <div>toleranser enligt:</div> | <div>Filnamn:</div> <div>Catalyst 5 med EasyWheel.icd</div>                         |                      | <div>Senast sparad:</div> <div>2017-12-06 14:41</div> |                                |
|                                                                                                                                | <div>Skala:</div> <div>1:2</div>                                                                                                                                        | <div>ISO 2768-1</div> <div>medium</div>            | <div>Vikt:</div> <div>2,129</div>                                                   | <div>Material:</div> | <div>Skapad datum:</div>                              |                                |
| <div> decon wheel</div>                     |                                                                                                                                                                         |                                                    | <div>Catalyst 5 / Catalyst 5Ti /</div> <div>Catalyst 5VX</div> <div>EasyWheel</div> |                      | <div>Bladnr:</div> <div>1 ( 1 )</div>                 | <div>Rev.</div> <div>...</div> |
|                                                                                                                                |                                                                                                                                                                         |                                                    |                                                                                     |                      | <div>Ritningsnr:</div> <div>VH3006-10</div>           |                                |