




Etac - Cross 5 & 6

MED3404




1	93440-120520	1	
2	XOC-25383-00	1	






1	HW069	1	
2	405242	1	
3	100 400 017	1	
4	MED2000	1	
5	100 400 004	1	
6	MED3100-01	1	
7	MED3100-04	2	
8	100 400 084	1	
9	MED3100-02	1	
10	XSM8X45	1	
11	100 400 001	1	
12	MED3100-03	2	
13	100 400 083	1	
14	XSM8ML	1	

Cross 5




Cross 6







NiMH




Li-ion



NiMH




Li-ion



NiMH




Li-ion



NiMH




Li-ion




NiMH



Li-ion



M16x1	30 Nm
M18x1	40 Nm



M3	1.2 Nm
M4	3 Nm
M5	6 Nm
M6	10 Nm
M8	25 Nm
M10	50 Nm
M12	85 Nm