




BeRollka - Findus

MED3171




1	93440-120520	1	
2	XOC-25383-00	1	





1	HW069	1	
2	405239	1	
3	MED2000	1	
4	100 400 017	1	
5	100 400 004	1	
6	100 400 036	2	
7	XSM16FJ	1	
8	XSM16M	1	
9	XSM8X40	1	
10	XSM8BR	2	
11	100 400 001	1	
12	100 400 037	2	
13	XSM8ML	1	



NiMH




Li-ion



NiMH




Li-ion



NiMH



Li-ion



Option


(Separate assembly instruction)



NiMH




Li-ion



NiMH




Li-ion




NiMH



Li-ion



M16x1	30 Nm
M18x1	40 Nm



M3	1.2 Nm
M4	3 Nm
M5	6 Nm
M6	10 Nm
M8	25 Nm
M10	50 Nm
M12	85 Nm