











Meyra - Vario 1.750
MED3096





1	MED2000	1	
2	XSM16BR	1	
3	100 400 017	1	
4	100 400 007	2	
5	HW069	1	
6	405239	1	
7	XSM16FJ	1	
8	XSM16M	1	
9	XSM8BR	1	
10	XSM8X40	1	
11	XSM8ML	1	


NiMH



Li-ion




Code 107
MED4053R



1	MED1590R	1	 
2	XSM6X16	2	
3	XSM6ML	2	
4	XOC-48160-B0	1	

 = MED4053L

Code 70
MED4037R



1	XSM6ML	2	
2	XSM6BR	2	
3	XSM6X14	2	
4	MED4036R	1	
5	X0C-48160-B0	1	



= MED4037L




Code 70
MED4038R

1	XSM6ML	2	
2	XSM6BR	2	
3	XSM6X14	2	
4	MED4036R	1	
5	X0F-48160-D00	1	



= MED4038L

Code 577
MED4071R




1	XSM6X14	2	
2	X0C-48160-B0	1	
3	MED1662	1	
4	XSM6BR	2	
5	XSM6ML	2	



= MED4071L




Code 577
MED4072R



1	XSM6X14	2	
2	X0F-48160-D00	1	
3	MED1662	1	
4	XSM6BR	2	
5	XSM6ML	2	



= MED4072L




1	XSM6BRTG	1	
2	XOD-2612K-00	1	
3	XA2-2612E-00	1	
4	XSM6BR	1	
5	MC6SM6X50	1	



MED1218 + MED1643



M3	1,2 Nm
M4	3 Nm
M5	6 Nm
M6	10 Nm
M8	25 Nm
M10	50 Nm
M12	85 Nm



M16x1	30 Nm
M18x1	40 Nm