


decon
mobility



ADVENTUS



Conquer any terrain with confidence – Adventus transforms your manual wheelchair into an off-road companion.

Its robust 10" front wheel lifts the small casters off the ground, preventing obstacles like roots, stones, or thresholds from stopping you. The result is a smoother, more stable, and safer ride on forest paths, gravel, cobblestones, sandy beaches, or uneven pavements.

Adventus is quick to attach and detach with its tool-free installation, and when not in use, it stores easily. Lightweight yet durable, it is compatible with most manual wheelchairs, both rigid and with foldable/folding frames, and gives users more freedom without compromising comfort or control.

Whether enjoying a walk in nature, joining family at the beach, or navigating uneven city streets, Adventus makes every journey easier. More than just an accessory, it opens up new opportunities to explore the world on your own terms.

Adventus – making the world more accessible, step by step.

Specifications

Length	ca. 450 mm
Height	ca. 350 mm
Wheel dimension	10"
Weight	2,45 kg
CE marking	In compliance with MDR 2017/745 class 1



Please scan the QR-code for more information.



Foldable frame chair adapter

The Foldable Frame Adapter gives your Adventus wheelchair maximum flexibility and compatibility. With this smart solution, Adventus can easily be mounted on wheelchairs with a cross frame design, ensuring a secure and stable fit. The adapter makes it possible to enjoy all the benefits of Adventus without compromising on safety or comfort, even if your wheelchair has a folding frame.

Whether for everyday use or more demanding environments, the Foldable Frame Adapter guarantees that Adventus integrates seamlessly, offering reliability and ease of use.



Drive easily over uneven ground

With Adventus, you can move confidently across uneven surfaces such as gravel, cobblestones, or grass. The strong and durable design improves traction and stability, reducing the effort required to push the wheelchair. This makes outdoor mobility smoother, safer, and more enjoyable – so you can focus on where you're going, not the obstacles in your way.



Please scan the QR-code for more information.

Together we can create freedom!

Follow us:



decon.se

